



Step-by-Step Communicator

What is a step-by-step?

This device allows for an individual to record any series of messages (up to 2 minutes) in to the Step-by-step communicator. You press the button surface one time to hear the first message. Hit the button again to hear the second message within the sequence, and so on. Within the 2 minutes of allotted time, you can record as many messages as needed.

This device has an angled base with a 2.5-inch activation surface. It requires a 9-volt battery that lasts for weeks to months depending on usage.

Instructions to record:

1. Select the desired level using the “level” switch located on the side of the device. Press and HOLD the “Record” button until the light activates and then release the button.
2. Press and hold the switch button (on the top) and begin speaking once you hear the beep. Note: The record light will flash while recording
3. Release when finished.
4. Repeat for additional messages within the sequence.
5. Press and release the record button to exit recording mode.
6. Your step-by step is now ready to use.

Ways to use the Step-by-Step:

- **About me:** My name is...., I’m _____ years old, My favorite song is....., etc.
- **Sing a song:** Itsy bitsy spider went up the water spout -press button for more--down came the rain..
- **Personal Narrative:** “Hi, daddy! Do you want to hear about my day?”
- **Games:** Simon says, Red-light Green-light
- **Ask and answer questions:** Who, What, Where, When, Why
- **Nursery rhymes/poems:** B-I-N-G-O
- **Counting:** Hide and Seek

For more information and pricing details: <https://www.ablenetinc.com>

Follow this QR Code to a YouTube video about using Step-by-Step Communicators

